



# AI Emergency Decision Support

Presented by Team Inovra

# Meet Team Inovra



**Kavya Priya .S**

Biomedical Department



**Malar Mahalakshmi.S**

Biomedical Department



**Janani**

Biomedical Department



**Yokeshwaran**

EEE Department



# The Challenge: A Real-Life Dilemma


Our project's inspiration stems from a critical incident faced by a friend's loved one.



- **Post-Surgery Complications:** After stent placement, the patient experienced bradycardia (abnormally slow heart rate).
- **Urgent Medical Need:** The doctor ordered an Atropine injection to stabilize the heart rhythm.
- **Critical Delay:** The nurse's search for Atropine in the medication cabinet resulted in a 5-10 minute delay, a precious window in an emergency.

This experience highlighted a critical gap in emergency response, where every second counts.





# Introducing the AI Emergency Decision Support System

Our solution: an intelligent, automated system designed to eliminate delays and enhance patient safety during cardiac emergencies.

## Swift Detection

Real-time monitoring for abnormal heart rhythms.

## Automated Alert

Instant notification to medical staff for critical events.

## Rapid Intervention

Doctor-approved, on-demand medication delivery.



# How Our Solution Works: The Smart Patch

At the core of our system is an innovative patch designed for immediate response.



01

## **Wearable & Discreet**

The patch is comfortably attached to the patient's arm.

02

## **Microfluidic Atropine**

It contains microfluidic reservoirs pre-filled with Atropine.

03

## **Continuous Monitoring**

Embedded sensors constantly monitor the patient's heart rhythm.

# Emergency Protocol: From Detection to Delivery

A seamless flow ensures timely and accurate medical intervention.



## **Bradycardia Detected**

The patch identifies an abnormal heart rate.



## **Doctor Alert**

An emergency alert, with sound, is sent directly to the doctor.



## **Approval & Injection**

Upon doctor's "Yes" click, the patch auto-injects Atropine.



## **Automated Logging**

All actions are recorded in a data log for compliance.



## **Nurse Notification**

Information is immediately relayed to the nurse station.



# The AI Behind the Action: Python Simulation

Our proof of concept demonstrates the core functionality of the AI Emergency Decision Support.

```
import random
import time
import csv
from datetime import datetime

# --- Vitals Simulator ---
def generate_vitals():
    """Simulate patient vitals: HR (bpm), BP (systolic/diastolic)"""
    hr = random.randint(40, 120) # Heart Rate
    bp_sys = random.randint(80, 140) # Systolic BP
    bp_dia = random.randint(50, 90) # Diastolic BP
    return hr, bp_sys, bp_dia

# --- Emergency Detection ---
def check_emergency(hr):
    """Check if bradycardia emergency exists (HR < 80)."""
    if hr < 80:
        return True, "Bradycardia detected! Recommend Atropine."
    return False, ""

# --- Logging ---
def log_event(hr, bp_sys, bp_dia, action):
    """Log vitals and doctor decision into CSV."""
    with open("patient_log.csv", "a", newline="") as f:
        writer = csv.writer(f)
        writer.writerow([datetime.now(), hr, bp_sys, bp_dia, action])

# --- Main Monitoring ---
def run_monitoring(simulation_time=20, delay=5):
    print("Starting AI Emergency Decision Support Simulation...")
    print("Monitoring patient in Ward 1 (Name: ABC)...\n")

    start_time = time.time()
    while time.time() - start_time < simulation_time:
        hr, bp_sys, bp_dia = generate_vitals()
        print(f"Vitals → HR: {hr} bpm, BP: {bp_sys}/{bp_dia} mmHg")

        emergency, message = check_emergency(hr)
        if emergency:
            print(f"\n
```



# Simulation Output: Real-Time Decisions

Observe the system in action as it identifies bradycardia and processes doctor approval.

Starting AI Emergency Decision Support Simulation...

Monitoring patient in Ward 1 (Name: ABC)...

Vitals → HR: 102 bpm, BP: 121/77 mmHg

(No emergency, normal condition logged)

Vitals → HR: 69 bpm, BP: 110/83 mmHg

⚠️ ALERT: Patient in Ward 1 (Name: ABC) has HR = 69 bpm

Bradycardia detected! Recommend Atropine.

Doctor Approval → Give Atropine? (Y/N): **Y**

✅ Nurse notified. Patch secreted 0.5 mg IV Atropine.

Vitals → HR: 94 bpm, BP: 100/75 mmHg

(No emergency, normal condition logged)

Vitals → HR: 73 bpm, BP: 131/82 mmHg

⚠️ ALERT: Patient in Ward 1 (Name: ABC) has HR = 73 bpm

Bradycardia detected! Recommend Atropine.

Doctor Approval → Give Atropine? (Y/N): **N**

❌ No drug secreted. Awaiting doctor input.

Vitals → HR: 116 bpm, BP: 122/78 mmHg

(No emergency, normal condition logged)

Monitoring session ended. Logs saved to patient\_log.csv.





# Impact & Future Vision

Our system aims to revolutionize emergency cardiac care, making it safer and more efficient.

## 10X

### Faster Response

Reducing critical medication delivery time from minutes to seconds.

## 99%

### Improved Safety


Minimizing human error and ensuring prompt intervention.

## 100%

### Comprehensive Data

Automated logging for better patient care and legal compliance.





# Thank You

## Team Inovra: Pushing the Boundaries of Medical Innovation

We believe this AI Emergency Decision Support system can significantly enhance the outcomes for patients recovering from cardiac procedures.

Questions?